

FarView



FarView

Portable Video Magnifier



Accelerating Innovations with the New Optelec FarView

Designed for the active lifestyle, Optelec's FarView offers the benefits of distance and near viewing in an ultra compact and stylish design. Using a center mounted auto-focus camera the FarView allows you to access



information with ease. Use the FarView on the road to read street names, departure times, or to select your favorite coffee. The FarView is also a great tool to use during a business meeting or at school to capture essential information on the whiteboard.

The FarView gives you the independence to read anywhere. You can capture and store your favorite book or magazine and read at your leisure. The auto scroll feature combined with intuitive functionality makes the FarView easy to use. With multiple color options and contrast modes the FarView can be customized to enhance any reading experience.

**To learn more about the FarView,
view videos and read additional information
visit us at www.Optelec.com.**



10. Up to 4 hours



Distance viewing



Reading



Multiple viewing options



USB & VGA connections



-2X-90X Zoom



10.2 ounces



6.6 x 3.2 x 1.3 in.



Memory storage
100 images



Image capture

Optelec U.S., Inc. Corporate Headquarters

Tel: 800-826-4200 • Fax: 800-368-4111

www.Optelec.com

