

MultiView



MultiView™

Portable Video Magnifier



Near, Distance and Self Viewing Solution

The Optelec MultiView is a transportable, multipurpose CCTV for document reading, distance viewing and self viewing. A versatile 3-in-1 camera and 19" widescreen display has a small footprint, allowing for simple transit and easy storage.

The high performance camera can be easily operated with the left or right hand. Additionally, the light weight camera can effortlessly be disconnected from the arm and moved to multiple work stations. Simply transport the MultiView from bedroom to home office or from home to school.

As a standard, the MultiView offers the Optelec smooth gliding XY reading table for ergonomic reading resulting in decreased eye fatigue. As an added benefit, the reading table is conveniently stored behind the unit when not in use.

Operation of the MultiView is simple with the intuitively placed control panel on the camera head. An advanced control panel comes standard with all units. For increased reading speed and comfort, the advanced control panel can be placed on either side of the unit, allowing the user to concentrate on the screen.

The Optelec MultiView is the perfect low vision solution for those with multiple needs. Read, write, self view and more with the Optelec MultiView!



Document reading



3.5-79x zoom



Self viewing



Distance viewing



Writing



19" widescreen



Multiple viewing options



55 x 48 x 21 cm
21.7 x 17.3 x 7.9 in



11.9 kg/26 lbs



Snapshot

Optelec U.S., Inc. Corporate Headquarters
Tel: 800-826-4200 • Fax: 800-368-4111
www.Optelec.com

Optelec
Tieman Group
Accelerating Innovation