


Generation 2



ClearView+[™]
Desktop Video Magnifier

Optelec
Life is worth enjoying

The Generation 2 ClearView+ sets the standard for desktop reading machines. Re-engineered for maximum comfort, functionality and image quality, nothing on the market compares.



Redesigned using state of the art HQV® image processing technology, the image quality is unparalleled. The all new ClearView+ has smoother character edges, vivid colors and higher contrast.



All units are HD compatible and can be upgraded to true HD image quality, enabling full A4/Letter text overview and very smooth characters, even at very low magnification levels. The advanced technology creates an optimal reading environment with smooth, crisp images, decreasing eye fatigue. This means faster reading for longer periods of time.



Already renowned for simplicity, the Generation 2 takes the concept one step further. The ClearView+ can be upgraded with an “AutoZoom” feature pack or the “Zero Button” module. These upgrades automatically convert and display text in the user’s preferred pre-set magnification and contrast. Simply rest your hands on the comfortable wrist pad, place your document on the xy-table and the ClearView+ takes care of the rest.



The ClearView+ redefines ergonomic design. For perfect ergonomic distance between the eye and the monitor, the ClearView+ 19” and 22” monitors come with a five directional “Ultra Flexible TFT Arm”. The Ultra Flexible Arm (UFA) allows the user to adjust the monitor to their preferred position with ease. The screen can be moved up or down, telescoped forward or backward, tilted, and angled left or right. The UFA also rotates 90-degrees and automatically flips the image to enable column reading.

No matter what you would like to do – view photos, read your favorite book or write letters – the Optelec ClearView+ can be fit to meet your needs.



Writing




Reading



Multiple Viewing



Lines & Windows



3X-65X Zoom



46.3-65 pounds



19.5 x 16.5 x 26.5 in.

Optelec U.S., Inc. Corporate Headquarters

Tel: 800-826-4200 • Fax: 800-368-4111

www.Optelec.com

Optelec
Life is worth enjoying