

**TV Glasses**



**Optelec**  
**Power TV Glasses**

[www.Optelec.com](http://www.Optelec.com)

# TV Glasses



## Optelec Power TV Glasses

Lightweight, adjustable Optelec Power TV Glasses provide 2.1X magnification for distance viewing. These comfortable glasses are perfect for watching television, movies, sporting events or other outdoor activities. Using a binocular telescope system, the Optelec Power TV glasses help to maximize the size of the television screen and other objects.

With adjustable lenses, the glasses can easily adjust and the hands free frames allows for complete enjoyment of multiple activities.



## Features

- 2.1X magnification
- 1.7 ounces in weight
- Twin lens design
- Hands-free
- 20 degree field of view
- Focal distance of 10 feet

Optelec US  
3030 Enterprise Court, Suite C  
Vista, California 92081  
800-826-4200  
[www.Optelec.com](http://www.Optelec.com)

**Optelec**  
Life is worth enjoying